roeko activ point

Instructions for use	EN
----------------------	----

Definition:
ROEKO *activ points* are time release preparations, which release chlorhexidine diacetate from a guttapercha matrix.

Composition:
activ points are made of chlorhexidine diacetate (approx. 5 %), guttapercha, ZnO, BaSO4 and colouring agents.

Indications:

- for temporary filling of root canals to prevent reinfection
- for emergency endodontic treatment

Properties:
The point is ready for use, firm yet flexible for easy introduction into root canals. Pure chlorhexidine diacetate is homogeneously distributed throughout the guttapercha matrix. The ISO points are ISO standard, dimensionally stable and orange to prevent confusion with gutta-percha or Calcium Hydroxide points.

Release of chlorhexidine:
Chlorhexidine diacetate is released into the fluid which flows into the canal through the dentine tubules and apex immediately after drying. The point allows sufficient chlorhexidine diacetate to be applied effectively into the canal.

Application:
An *activ point* corresponding to the last used root canal instrument, or one size smaller, should be used. The predetermined length should be marked on the point before applying it down to the apex with tweezers. The point should sit passively in the canal without any condensation. In oval or extremely conical canals extra points can be applied. Usually the access opening is large enough to fold in the excess length. This also facilitates removal later. All standard materials can be used to close the opening.

Duration:
The guttapercha matrix allows large quantities of chlorhexidine diacetate to be released from the surface of the point. Chlorhexidine diacetate which is embedded in the point is released more slowly. Its affinity to dentine means it remains present for some time. It is recommended that the point is left in the canal for 1-3 weeks. Afterwards the point should be replaced or the canal be permanently filled. In certain clinical situations the point should be replaced at shorter intervals (after 2-3 days).

Removal of the ROEKO *activ point*:
The release of chlorhexidine diacetate into the surroundings does not affect the stability of the point. It remains stable over a longer period of time and is easy to remove with tweezers or a probe.

Side effects:
Contact with mucous membranes may result in short term irritation and /or a reduced sense of taste. A reversible brown discoloration of the root canal may occur.

Contraindications:
In case of a known chlorhexidine diacetate allergy *activ point* should not be used.

Interactions:
If sodium hypochlorite or alcohol are used to rinse the canal during preparation there may be a synergetic effect on using chlorhexidine i.e. the efficacy may be increased, but there is no risk for the patient. Therefore the points may be used in a canal which has not been fully dried.

Use and storage:

- for single use only
- store in a dry, cool place (not above 25 °C)
- do not autoclave, do not sterilize in dry heat
- shelf life: 3 years
- this medical device should not be used after the expiration date

Warning:
Medical devices should be kept out of reach of children.

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Coltène/Whaledent GmbH + Co. KG